




International Filter Products

Revolutionary systems for your industry

Nominal Dimension: All dimensions are in Millimeter

Part Number: LI1DS065S020G4H4H-ETO



Hold-Up Volume (Approx.):

Upstream Volume: 916 ml

Downstream Volume: 127 ml

Typical Filtered Volume:

100.0L – 400.0L

Sterilization:

Capsules can be autoclaved 25 times at 125 °C (257°F) for 30 minutes or chemically sanitized in situ using common sanitizing agents or hot water at 90 °C (194°F) for a limited time (dependent on time and temperature).

Warning: The filters cannot be sterilized by steam-in-place (SIP)