




International Filter Products

Revolutionary systems for your industry

Nominal Dimension: All dimensions are in Millimeter

Part Number: LZT3XS020GTCTC



Hold-Up Volume (Approx.):

Upstream Volume: 2.28 L

Downstream Volume: 558.40 mL

Typical Filtered Volume:

300.0L – 1,200.0L

Filter Integrity:

The finished product was sampled and shown to exhibit a minimum bubble point of ≥ 18 psi (1.2 bar) in 60%IPA / 40% water and ≥ 50 psi (3.5 bar) in DI Water (at 22°C).

Sterilization:

Gamma Irradiation: The filters can be Gamma Sterilized up to 50kGy.

Autoclave: The filters can be autoclaved up to 10 cycles at 125°C (257°F) for 30 minutes.

Warning: The filters cannot be sterilized by steam-in-place (SIP)

Shelf Life:

The LZ capsules have a shelf life of 3 years from the date of manufacture.



7914 Ajay Drive, Sun Valley, CA 91352
Email: sales@internationalfilterproducts.com
Phone: 818-504-8115
www.internationalfilterproducts.com