




Sol-M[®] Standard Insulin Syringe – U100

- ▶ The Sol-M[®] Standard Insulin Syringe with fixed needle is used to inject insulin into the body. The syringe has a low dead space design which results in less medication waste and improved dosing accuracy.



Features

- 1 Tri-bevel, anti-coring and silicone oil lubricated**
Makes injection more comfortable
- 2 Low dead space**
Results in 80% less medication waste compared to the ISO standard for conventional insulin syringes
- 3 Transparent barrel with bold scale**
Easy to read units for dosing accuracy
- 4 Large flange**
Offers secure grip and easier medication dispensing



Sol-M® Standard Insulin Syringe - U100

Sol-M® Standard Insulin Syringe - (PE Bag - 10 pcs)

Ref. Number	Volume	Gauge	Needle Length (mm)	Needle Length (in)	Box Qty	Case Qty
1612712B	1mL	27G	12.5mm	1/2"	100	1,000
1612758B	1mL	27G	16mm	5/8"	100	1,000
1612812B	1mL	28G	12.5mm	1/2"	100	1,000
1612912B	1mL	29G	12.5mm	1/2"	100	1,000
1613012B	1mL	30G	12.5mm	1/2"	100	1,000
1613051B	1mL	30G	8mm	5/16"	100	1,000
1613151B	1mL	31G	8mm	5/16"	100	1,000
161311564B	1mL	31G	6mm	15/64"	100	1,000
1652712B	0.5mL	27G	12.5mm	1/2"	100	1,000
1652812B	0.5mL	28G	12.5mm	1/2"	100	1,000
1652912B	0.5mL	29G	12.5mm	1/2"	100	1,000
1653012B	0.5mL	30G	12.5mm	1/2"	100	1,000
1653051B	0.5mL	30G	8mm	5/16"	100	1,000
165311564B	0.5mL	31G	6mm	15/64"	100	1,000
1653151B	0.5mL	31G	8mm	5/16"	100	1,000
1633012B	0.3mL	30G	12.5mm	1/2"	100	1,000
1633051B	0.3mL	30G	8mm	5/16"	100	1,000
163311564B	0.3mL	31G	6mm	15/64"	100	1,000
1633151B	0.3mL	31G	8mm	5/16"	100	1,000

For more information, please contact:



Sol-Millennium Medical Inc.
311 S Wacker Dr. Suite 4100
Chicago, IL 60606



833-637-0857



info@solm.com

NACA-DC-463-0425