

Sol - M™ Hypodermic Needle

- Sol-Millennium® offers a wide range of hypodermic needles with an emphasis on ease of use and patient comfort.

Features:

- Triple-bevel lubricated stainless steel bevel engineered to minimize patient discomfort
- Color-coded hubs for facilitating gauge identification
- Compatible with all standard Luer Lock and Luer Slip Tip syringes (ISO 80369-7)
- Designed according to ISO 7864
- Not manufactured with Latex, PVC and DEHP



CE 2797

Color-coded and translucent hubs

Sol-M™ Hypodermic Needle

Ref. Number	Product Description	Gauge	Length (inches)	Length (mm)	Box Qty	Case Qty
111615	Sol-M Hypodermic Needle 16G x 1 1/2"	16G	1 1/2"	40 mm	100	1000
111810	Sol-M Hypodermic Needle 18G x 1"	18G	1"	25 mm	100	1000
111815	Sol-M Hypodermic Needle 18G x 1 1/2"	18G	1 1/2"	40 mm	100	1000
111820	Sol-M Hypodermic Needle 18G x 2"	18G	2"	50mm	100	1000
111910	Sol-M Hypodermic Needle 19G x 1"	19 G	1"	25 mm	100	1000
111915	Sol-M Hypodermic Needle 19G x 1 1/2"	19 G	1 1/2"	40 mm	100	1000
111920	Sol-M Hypodermic Needle 19G x 2"	19 G	2"	50 mm	100	1000
112010	Sol-M Hypodermic Needle 20G x 1"	20G	1"	25 mm	100	1000
112015	Sol-M Hypodermic Needle 20G x 1 1/2"	20G	1 1/2"	40 mm	100	1000
112110	Sol-M Hypodermic Needle 21G x 1"	21G	1"	25 mm	100	1000
112115	Sol-M Hypodermic Needle 21G x 1 1/2"	21G	1 1/2"	40 mm	100	1000
112120	Sol-M Hypodermic Needle 21G x 2"	21G	2"	50 mm	100	1000
112210	Sol-M Hypodermic Needle 22G x 1"	22G	1"	25 mm	100	1000
1122125	Sol-M Hypodermic Needle 22G x 1 1/4"	22G	1 1/4"	32 mm	100	1000
112215	Sol-M Hypodermic Needle 22G x 1 1/2"	22G	1 1/2"	40 mm	100	1000
112310	Sol-M Hypodermic Needle 23G x 1"	23G	1"	25 mm	100	1000
1123125	Sol-M Hypodermic Needle 23G x 1 1/4"	23G	1 1/4"	32 mm	100	1000
112410	Sol-M Hypodermic Needle 24G x 1"	24G	1"	25 mm	100	1000
112506	Sol-M Hypodermic Needle 25G x 5/8"	25G	5/8"	16 mm	100	1000
112510	Sol-M Hypodermic Needle 25G x 1"	25G	1"	25 mm	100	1000
112612	Sol-M Hypodermic Needle 26G x 1/2"	26G	1/2"	12.5 mm	100	1000
112658	Sol-M Hypodermic Needle 26G x 5/8"	26G	5/8"	16 mm	100	1000
112705	Sol-M Hypodermic Needle 27G x 1/2"	27G	1/2"	12.5 mm	100	1000
112734	Sol-M Hypodermic Needle 27G x 3/4"	27G	3/4"	20 mm	100	1000
112905	Sol-M Hypodermic Needle 29G x 1/2"	29G	1/2"	12.5 mm	100	1000
113005	Sol-M Hypodermic Needle 30G x 1/2"	30G	1/2"	12.5 mm	100	1000

Please contact your local SOL-MILLENNIUM® representative for additional sizes of the Sol-M™ Hypodermic Needle.

Intended use: Sol-M™ Hypodermic Needle is used to withdraw blood or other fluids from the body as well as inject medications into the body.

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/ injury.
- Sharp object. Risk of sharp injuries.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMA-MD-453-0423

